



Quick User's Guide – Lifeina Breathe Oximeter

This pulse oximeter is a non-invasive device intended for spot-checking oxygen saturation of arterial hemoglobin (SpO₂) and Pulse Rate of adults, adolescents and child patients in hospitals, hospital-type facilities and homecare.

1. Install two AAA batteries by opening compartment under device and following the +/- signs in the compartment.
2. Place one of your fingers into the rubber opening of the pulse oximeter. The finger most used is the index finger. Please ensure that you use the same finger every time to get comparable results.
3. Press the O/I button on the front panel to turn the pulse oximeter on. By repeatedly pressing the O/I button, you can turn the display or change display mode.
4. Keep your hands still during measurement, which should last about 5 seconds.
5. The result is displayed on the front screen.
6. Pressing the power switch for longer than one second will adjust the brightness of the oximeter. There are 10 levels of brightness. The default is level four.



Ideal position for your oximeter

Note: If the screen displays “?” the signal is unstable. Please keep your hands still and retry.

Discover more incredible medical devices on
www.lifeina.com