

# LifeinaBox



*For the optimum storage  
of fragile medication*

## **THE FREEDOM TO TRAVEL**

Designed for the safe storage of fragile medication at an optimum temperature (2 - 8°C). Keeps your medication cool and at a constant temperature at work, at home or in your car.

## **CONSTANT TEMPERATURE**

Refrigerates up to 8 standard injectable medicine pens or vials\*. Includes both a universal 110-240V power cord and a 12V DC power adapter for your car.

## **BETTER CONTROL**

The Lifeina app will monitor the temperature of your medication in real time and will send you automatic reminders when it is time for you to take your medication.

## **BUILT-IN BATTERY**

A 6-hour battery is included inside LifeinaBox, and separate plug-in battery packs are available.

## **RECOMMENDED USES**

Insulin, growth hormones, adrenalin, monoclonal antibodies, vaccines, interferons\*

[www.lifeina.com](http://www.lifeina.com)

\*For a comprehensive list of compatible pens and medications, log-on to [www.lifeina.com](http://www.lifeina.com)



*"It's time to start living the life you've imagined"*

Henry James

- ▶ *Compact and portable*
- ▶ *Connected to the Lifeina app for real-time monitoring of temperature*
- ▶ *Automatic medication reminders*
- ▶ *Low-noise fan for quiet operation*
- ▶ *Built-in battery*
- ▶ *Useable in car (12V) or home (110-240V)*

Did you know? The LifeinaBox mini-refrigerator uses thermoelectric technology. Thermoelectrics convert thermal energy into electrical energy to produce cold.

Thermoelectric refrigeration uses a principle called the "Peltier" effect to pump heat electronically. The Peltier effect is named after Jean Charles Athanase Peltier, a French scientist who discovered it in 1834.

*Live your life*



COMPACT DESIGN



CONNECTED



TRACK TEMPERATURES



6-HOUR BATTERY



USEABLE ANYWHERE

## Specifications

|                            |   |
|----------------------------|---|
| Dimensions:                | 220 x 107 x 110mm   |
| Cavity size:               | 180mm wide x 65mm deep x 40mm high (468cm <sup>2</sup> )                |
| Capacity:                  | 8 prefilled insulin pens  |
| Net Weight:                | 1200gm  |
| Voltage:                   | 12V DC auto socket cord or 110-240V AC domestic power cord, UL approved |
| Battery                    | 2 x 3.7V Li-ion 5200mAh   |
| Average power consumption: | 25 watts max  |
| Temperature range:         | 2°C to 8°C, +/- 1°C.  |
| Insulation:                | CFC-free polyurethane foam  |
| System:                    | Wear-free thermoelectric technology with interior semiconductor fan     |
| Material:                  | ABS and brushed aluminum  |
| Reference:                 | LifeinaBox  |

*Specifications are subject to technical improvement, changes or modifications without notice*



**Lifeina SAS**  
6 Rue de Castellane  
75008 Paris, France  
[hello@lifeina.com](mailto:hello@lifeina.com)  
[www.lifeina.com](http://www.lifeina.com)